



NE1D-2/2/2



Triplebath™ features three independently controlled unstirred water baths in one unit. Typically suitable for classroom use performing PCR routines in shallow temperature controlled water, each tank can be set to a different temperature for denaturing, annealing and elongation reducing solvent and consumables.

Comprising of high quality corrosion resistant stainless steel tanks, with a bright clean finish housed in a durable paint finished outer case. Painted surfaces feature anti-bacterial finish with good chemical resistance. 4 Litre chambers are fitted with side lifting handles for ease of transportation.

Each chamber comes complete with a removable “low” height stainless steel perforated shelf.

NE1-2/2/2      NE1-4/4/4      NE1D-2/2/2      NE1D-4/4/4

GENERAL				
Chamber Capacity	3 x 2 Litre	3 x 4 Litre	3 x 2 Litre	3 x 4 Litre
Temperature Range	Ambient + 5°C - 99°C			
Temperature Control	Analogue dial pointer, graduated in 5°C divisions		Digital, LED display showing actual or set temperature, resolution 1°C	
Sensitivity/Uniformity	±0.5°C / ±0.5°C		±0.25°C / ±0.1°C	
Heater Power (220-240V)	3 x 200W	3 x 400W	3 x 200W	3 x 400W
Voltage	230V / 50-60Hz			
DIMENSIONS				
External (mm)	535w x 332d x 135h	535w x 332d x 290h	535w x 332d x 135h	535w x 332d x 290h
Internal (mm)	300w x 150d x 65h	300w x 150d x 150h	300w x 150d x 65h	300w x 150d x 150h
Water Level Max. (mm)	30	130	30	130
Water Level Min. (mm)	10	55	10	55

# TripleBaths™

## Accessories

Please note : All lids can be interchanged.



### Stainless Steel Lift Off Gable Lid

Bath	NE1-2/2/2	NE1-4/4/4	NE1D-2/2/2	NE1D-4/4/4
Accessory Code	SL1-4 / SL1-4 / SL1-4			



### Stainless Steel Flat Lid

Bath	NE1-2/2/2	NE1-4/4/4	NE1D-2/2/2	NE1D-4/4/4
Accessory Code	LD-4 / LD-4 / LD-4			



### Clear Lift Off Gable Lid

Bath	NE1-2/2/2	NE1-4/4/4	NE1D-2/2/2	NE1D-4/4/4
Accessory Code	GL1-4AF / GL1-4AF / GL1-4AF			



### Polypropylene Spheres (Pack of 200)

Bath	NE1-2/2/2	NE1-4/4/4	NE1D-2/2/2	NE1D-4/4/4
No. of Packs	1 Pack per chamber			
Accessory Code	BP0368			



### Stainless Steel Test Tube Racks

Bath	NE1-2/2/2	NE1-4/4/4	NE1D-2/2/2	NE1D-4/4/4
No. of Racks	N/A	1	N/A	1
Accessory Code	6870 : 26 holes x 17mm DIA			
	6871 : 16 holes x 26mm DIA			
	6872 : 36 holes x 13mm DIA			
	6873 : 18 holes x 19mm DIA			
	6900 : 12 holes x 32mm DIA			

When using the racks on a 2 litre chamber, remove 1 shelf from the rack.



### Bath Armor™

Bath	NE1-2/2/2	NE1-4/4/4	NE1D-2/2/2	NE1D-4/4/4
No. of packs	4 Litre - 1 x BA-5		N/A	
	4 Litre - 1 x BA-5			
	4 Litre - 1 x BA-5			
Accessory Code	BA-5 : 4 Litre Pack BA-12 : 12 Litre Pack		N/A	

## Universal Accessories

CAL	Calibration Record
BX0616	Draining Syphon