NE1B Analogue Boiling Bath Boiling Baths

Clifton Range®



Designed for applications that require continuous boiling. This bath also has the flexibility where the control can be regulated from a gentle simmer to a vigorous boil.

Comprising of a stainless steel tank housed in a durable paint finished outer case. Painted surfaces feature a tough anti-bacterial finish with good chemical resistance. Complete with a removable stainless steel perforated false base, drain tap and clear lift off lid. The bath, which is simmer-statically controlled, is ideally suited to prolonged periods of boiling.

Fitted with a sheathed immersion element. With the heater directly immersed in water it provides rapid boiling. A standard feature is a constant level device, which maintains a fixed level of liquid thereby eliminating the need to regularly top up the bath. Not fitted on NE1B-18W.

	NE1B-14	NE18-18W	NE1B-28
GENERAL			
Chamber Capacity	14 Litres	18 Litres	28 Litres
Temperature Range	Simmer to boiling at 100°C		
Temperature Control	Featuring an energy regulator whi	ch controls the amount of boiling fro	m a gentle simmer to vigorous boil
Heater Power (220-240V)		2500W	
Voltage	230V / 50-60Hz		
Water Level Alarm	Features low liqui	d level float device, heating autom	atically suspended
Sterilising Tray	×		×
Clear Lift Off Lid Supplied			·
Drain Tap Fitted		\checkmark	
Constant Level Device	\checkmark	×	\checkmark
DIMENSIONS		·	·
External (mm)	362w x 33	32d x 290h	538w x 332d x 290h
Internal (mm)	325w x 300d x 150h	325w x 300d x 200h	500w x 300d x 200h



Low liquid float device

A low liquid float device protects the heaters and a red warning light on the front panel alerts users to refill the bath.



NE1B-18W

Raise and lower sterilising tray along with pair of guide rail assemblies comes as standard. Provides a safe method of lifting sterilised items from boiling water. Excess hot water drains away.

Nickel-Electro Ltd., Oldmixon Crescent, Weston-super-Mare, North Somerset, BS24 9BL, United Kingdom t +44 (0)1934 626691 w www.nickel-electro.co.uk e clifton@nickel-electro.co.uk



NE1B Analogue Boiling Bath Accessories



Stainless Steel Lift Off Gable Lid

	Bath	NE1B-14	NE1B-18W	NE1B-28
5	Accessory Code	N/A	N/A	SL1-22SS

Stainless Steel Flat Lid

Bath	NE1B-14	NE1B-18W	NE1B-28
Accessory Code	LD-14	N/A	LD-22



Stainless Steel Hinged Gable Lid

Bath	NE1B-14	NE1B-18W	NE1B-28
Accessory Code	SL1-14	SL1-14H/SS	

Clear Hinged Gable Lid

	Bath	NE1B-14	NE1B-18W	NE1B-28
>	Accessory Code	Code GL1-14		GL1-22



22

clift

Stainless Steel Concentric Ringed Lid

	Bath	NE1B-14	NE1B-18W	NE1B-28
) [Accessory Code	SLR1-14	N/A	SI B1-22
	7.0000001y 000e	SLR2-14	11/7	OLITI-ZZ

SLR1-14: 4 x 114mm holes

SLR2-14: 6 x 83mm holes

SLR1-22: 6 x 114mm holes

Polypropylene Spheres (Pack of 200)

	Bath	NE1B-14	NE1B-18W	NE1B-28
	No. of Packs	2 Packs		4 Packs
	Accessory Code	BP0368		

NE1B Analogue Boiling Bath Accessories



Stainless Steel Test Tube Racks

	Bath	NE1B-14	NE1B-18W	NE1B-28
>	No. of Racks	3		5
		6	870 : 26 holes x 17mm DI	A
	Accessory Code	6	871 : 16 holes x 26mm DI	A
		6	872 : 36 holes x 13mm DI	A
			873 : 18 holes x 19mm Dl	A
		6	900 : 12 holes x 32mm DI	A

Universal Accessories

CAL

Calibration Record